Chirurgie *Surgery*

Quality Products
Made in Germany

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Dr. Ivo Agabiti, Pesaro, Italien
www.sonosurgery.it

Dr. Martin Dürholt,
Bad Salzuflen

Dr. Ernst Fuchs Schaller,
Thalwil, Schweiz

Dr. Stefan Neumeyer,
Eschlkam

Dr. med. dent. Marcus Striegel,
Dr. Thomas Schwenk,
Nürnberg

Prof. Dr. Klaus-Ulrich Benner,
Germering

Dr. Dr. Karl-Heinz Heuckmann,
Chieming

Scientific advice:

*Dr. Ivo Agabiti, Pesaro, Italy
www.sonosurgery.it*

*Dr. Martin Dürholt,
Bad Salzuflen, Germany*

*Dr. Ernst Fuchs Schaller,
Thalwil, Switzerland*

*Dr. Stefan Neumeyer,
Eschlkam, Germany*

*Dr. med. dent. Marcus Striegel,
Dr. Thomas Schwenk,
Nuremberg, Germany*

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Stand: Februar 2015

Some of the products and designations mentioned in the text are trademarked, patented or copyrighted.

The absence of a special reference or the sign ® should not be interpreted as the absence of legal protection.

This publication is copyrighted. All rights, also with regard to translation, reprint and reproduction (also in the form of extracts) are reserved. No part of this publication may be reproduced or processed using electronic systems in any form or by any means (photocopying, microfilm or other methods) without the written permission of the editor.

Colours and products subject to alterations. Printing errors excepted.

As at February 2015

	<p>Chirurgische Knochenfräser <i>Surgical bone cutters</i></p>		<p>Implantologie <i>Implantology</i></p>
4 – 14	<p>Knochenfräser aus Hartmetall <i>Bone cutters made of tungsten carbide</i></p>	64 – 66	<p>Universelle Pilotbohrer <i>Universal pilot burs</i></p>
15	<p>Knochenfräser aus Keramik <i>Bone cutters made of ceramics</i></p>		
16 – 17	<p>Knochenfräser diamantiert <i>Bone cutters diamond-coated</i></p>	67 – 71	<p>Zubehör <i>Auxiliaries</i></p> <p>Instrumentenständer <i>Bur blocks</i></p>
18	<p>Innengekühlte Instrumente <i>Instruments with internal cooling</i></p>	72 – 75	<p>Inserttrays und Sterilcontainer <i>Insert trays and sterilization container</i></p>
19 – 20	<p>Trepanbohrer aus RF-Stahl <i>Trepan burs made of stainless steel</i></p>	76	<p>Zubehör Schallspitzen <i>Auxiliaries for sonic tips</i></p>
	<p>Chirurgische Schallspitzen <i>Surgical sonic tips</i></p>	77	<p>Zubehör <i>Auxiliaries</i></p>
21 – 26	<p>SonicLine Chirurgie <i>SonicLine Surgery</i></p>	78 – 79	<p>Klinische Fälle <i>Clinical cases</i></p>
27 – 29	<p>Schallspitzen für die chirurgische Kronenverlängerung <i>Sonic tips for surgical crown extension</i></p>		
30 – 32	<p>Hubfeilen <i>Files for reciprocating handpiece</i></p>		
33 – 35	<p>TMC Extrusion <i>TMC Extrusion</i></p>		
	<p>MaxilloPrep</p>		
36 – 41	<p>MaxilloPrep Spread-Condense <i>MaxilloPrep Spread-Condense</i></p>		
42 – 51	<p>Angle Modulation System <i>Angle Modulation System</i></p>		
52 – 55	<p>Pipe Augmentation System</p>		
56 – 60	<p>MaxilloPrep Bone <i>MaxilloPrep Bone</i></p>		
61 – 63	<p>MaxilloPrep Bonefix <i>MaxilloPrep Bonefix</i></p>		

Chirurgie

Surgery

Komet® offers a wide variety of bone cutters made of tungsten carbide, ceramics, diamond instruments and "Miniflex" discs for bone treatments. The instrument shapes were designed to suit the requirements of different indications. State-of-the-art blade geometries ensure excellent cutting efficiency and smooth operation. The high-quality instruments made by Komet are synonymous for optimal clinical results and minimally invasive, efficient preparations. Special trepan burs and pilot burs have been designed for implantological treatments. Even in small diameters of just 1 mm, they help determining the direction and depth of an implant site, thus ensuring an efficient preparation. We are particularly proud of our most recently developed innovations: our bone cutters H254E, H162SL and H162ST.

Range of instruments:

- Bone cutters
- Miniflex bone saw
- Instruments with internal cooling
- Trepan burs
- Bone chip extractors
- Universal pilot burs
- Bur blocks for surgical instruments

Systems for bone extraction and augmentation

Komet offers special systems for oral surgery and pre-implantology, such as the MaxilloPrep Spread-Condense screws or the Angle Modulation system for minimally invasive bone spreading.

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


Bone cutter round, high-efficiency cutting design

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
Bone cutter round, high-efficiency cutting design, ZrN coated

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








Bone cutter with high-efficiency cutting design

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Bone cutter, cross cut tapered fissure long

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
Bone cutter, round end cross cut tapered fissure

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






Bone cutter, Lindemann, with special staggered toothing

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Bone cutter, Lindemann, with special staggered toothing, ZrN coated

Knochenfräser H162ST

Bone cutter H162ST

A new cut, a new generation: Sharp as the teeth of the proverbial sabre tooth tiger, Komet®'s new ST tothing makes work with the H162ST a completely new experience. This tapered bone cutter is particularly effective and cuts with unrivalled, perfect sharpness.

The advantages at a glance:

- Perfect sharpness
- Ideal cutting behavior
- Maximum control

The instrument is suitable for all sectors of dental surgery: From bone cuts during an osteotomy via hemisections and axial bone perforations to the surgical removal of retained teeth or root fragments.

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

Bone cutter tapered, maximum cutting performance thanks to special ST-tothing

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






*Bone cutter, Lindemann, with special staggered tothing,
ZrN coated*

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
*End-cutting bur for reducing bone substance during
surgical crown extension, for recreating the natural
biological width or for creating a flat preparation floor in
the cavity, laser marking at 4 mm*

K 160 A

		5	5	5
Größe · Size	Ø 1/10 mm	023	027	031


Winkelstück lang · Right-angle long (RAL)

K160A.205. ...

023 027 031

Handstück · Handpiece (HP)

K160A.104. ...

023 027 031

⊖_{max} 40000 min⁻¹/rpm

Gebrauchsmuster, Patente / Utility model, patents

EP 1 539 018*

* angemeldet / * pending

Knochenfräser rund, Keramik

Bone cutter round, made of ceramics

K 157

		5
Größe · Size	Ø 1/10 mm	016
L	mm	9,0


FG · Friction Grip (FG)

K157.314. ...

016

Winkelstück · Right-angle (RA)

K157.204. ...

016


Winkelstück lang · Right-angle long (RAL)

K157.205. ...

016

Handstück · Handpiece (HP)

K157.104. ...

016

⊖ = ⊖_{max} 40000 min⁻¹/rpm

■ = ⊖_{max} 160000 min⁻¹/rpm

Knochenfräser, Keramik

Bone cutter, made of ceramics

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Diamond-coated bone cutters

Contrary to cutting instruments, diamond coated cutters work in abrasive mode. The round shape can for example be used for mobilising a lateral bone window as part of a sinus lift.

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
Diamond disc Miniflex

The extra fine Miniflex diamond disc is most frequently used for splitting a jaw section to widen the alveolar ridge prior to a restoration with an implant. It is equally suitable for removing bone blocks in the region of the jaw angle. Use with disc guard.

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

Diamond cutter for apicectomy and separation of wisdom teeth

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
Miniflex diamond disc for bone-lid method

Apicectomy in the molar area, osteoplastic surgery of the maxillary sinus

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


Trepanbohrer, RF-Stahl

Trepan Burs, Stainless steel

The following three trepan burs with different designs are intended for the creation of different kinds of hollow perforations.

227A | Safe explantation

- High-efficiency toothing specially designed for cutting bone, laser depth markings to permit precise exposure of the implants to the required depth, large openings in the working part for a better view of the surgical site.

227B | Preparation of bone cylinders

- High-efficiency toothing specially designed for preparing bone cylinders with a diameter of 2.9 to 8.9, available in different lengths.

227C | Apical ectomy

- Suggested by Dr. Dr. Ayad, Münster, Germany
- Reduced trepanation depth of merely 5 mm.

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



Trepan bur for preparation of bone cylinders, stainless steel

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
Trepan bur for preparation of bone cylinders in the course of an apicectomy, stainless steel

Dr. Ivo Agabiti,
Pesaro, Italien
www.sonosurgery.it

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Sonic tips for surgery

Extremely fine incisions in bones

Developed in close cooperation with Dr. Ivo Agabiti from Pesaro, Italy, these sonic tips are a great leap forward in the quality of oral surgery. The innovative tips allow particularly gentle, conservative work on bones within the scope of restorations involving implants.

The sonic tips SFS100, SFS101 and SFS 102 are indicated for splitting of the crest and lateral incisions in the bone as part of bone spreading treatments and for gentle extraction of a tooth from its alveolar compartment.

The saucer-shaped sonic tips SFS103 and SFS104 as well as the oval SFS105 are particularly suited for gently detaching the sinus membrane within the course of an external sinus lift operation.

SFS 100/101/102:

- The extremely fine cuts of merely 0.25 mm guarantee maximum conservation of substance in the bone structure.
- The blade length of 10.7 mm provides sufficient penetration depth to allow vertical countersinking into the spongy bone; the width of the blade is 3.5 mm.
- Thanks to their oscillating movements, the tips are particularly gentle on the soft tissue.
- The sonic tips are used in oscillating mode in a suitable sonic hand piece, e.g. Komet® SF1LM, which is placed on a dental turbine. A special motor is not required. The operating frequency of the oscillating, elliptic movements is approx. 6,000 Hz.
- Sterile cooling solution is supplied via a cooling adapter equipped with small tubes that are screwed in between the sonic hand piece and the sonic tip.

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

Rostfreier Stahl

*Diamond coated, round, straight, for the preparation of a lateral window/
external sinus lift*

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Rostfreier Stahl

*Diamond coated, round, straight, fine grit, for the preparation of a lateral
window/external sinus lift*

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Dr. Marcus Striegel

Sonic tips SFS120/121/122

Minimally invasive surgical crown extension

The sonic tips SFS120, SFS121 and SFS122 facilitate minimally invasive surgical crown extension without the creation of a flap. These sonic tips were developed in close cooperation with Dr. Thomas Schwenk and Dr. Marcus Striegel, Nuremberg, Germany.

They carry out a surgical crown extension for example in cases where the biological width has been violated or when an asymmetric dental arcade or a "gummy smile" are to be rectified.

- *The tip is available in two diameters, allowing an optimal aesthetic treatment of the eye-catching areas. SFS120: Size 020 for teeth 2 and 4, size 030 for teeth 1 and 3.*

- *Posterior use: Even in the posterior region, the bone is often too close to the crown margin.*

- *The treatment is carried out in the tried and tested clinical manner: Having measured the patient's individual biological width, the treatment is planned with the help of a mock-up, wax-up or template. After the subsequent gingival correction, the osteotomy to recreate the biological width without creation of a flap is carried out with the sonic tip by slightly pivoting the tip in mesial and distal direction.*

- *To be used in a suitable sonic hand piece, e.g. Komet® SF1LM*

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Reciprocating files as part of the TissueMaster Concept®

Instruments performing stroke movements tend to be used in special treatments in the dental surgery, and the popularity of these treatments just keeps on growing.

Oscillating reciprocal files are frequently used for trimming restored interproximal surfaces, in orthodontics (e.g. stripping), in prophylaxis or in the preparation of root surfaces.

The innovative new reciprocal files developed by Komet give oscillating movements a whole new momentum. Developed in close cooperation with Dr. Stefan Neumeyer, these state-of-the-art files are ideally suitable for work on hard and soft tissue within the scope of Dr. Neumeyer's TissueMaster Concept, short TMC. Don't just take our word for it. Give these new files a try and see for yourself!

The components of the innovative range of reciprocal files:

- Diamond coated files for mechanical periodontal treatments
- Saw blades for cutting bone
- Scalpels for work on soft tissue - initial, advanced and extended versions with a perfectly ergonomic, curved shape

Oscillating, straight and effective - These innovative files mark the dawn of a new era of intuitive, tactile work right at your fingertips.

Advantages:

- Great results in little time
- Absolute precision
- Effective performance
- Outstanding intuitive control

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*Reciprocating file, periodontal, diamond coated, „Surfer“,
with applicator/support*

For use in the reciprocating handpiece, stainless steel

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*Reciprocating file, periodontal, diamond coated,
„Smoother“, with applicator/support*

For use in the reciprocating handpiece, stainless steel

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*Reciprocating file, periodontal, diamond coated, „Beaver“,
with applicator/support*

For use in the reciprocating handpiece, stainless steel

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Reciprocating file saw, „jigsaw“, with applicator/support

For use in the reciprocating handpiece, stainless steel

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*Reciprocating file saw, „hedgehog“, with applicator/
support*

For use in the reciprocating handpiece, stainless steel

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
Extrusion therapy with the TissueMaster Concept®

The extraction of a tooth launches resorption processes that lead to a considerable loss of periodontal and alveolar hard and soft tissue. The progression of these processes and the resulting tissue damage require therapeutic measures.

The orthodontic extrusion within the scope of the TissueMaster concept developed by Dr. Stefan Neumeyer constitutes a minimally invasive therapy with decisive advantages: A coronal movement of the adjacent periodontal and alveolar tissue is generated, which prevents a recessive loss of tissue.

The extrusion therapy is a useful complementary treatment with a view to optimising the creation of the implant bed. The basic steps are: "Extraction, replanting, extrusion, implanting".

Advantages:

- Replantation in order to maintain alveolar structures
- The aim of the extrusion therapy is to launch a coronal movement of the alveolar structures
- Logical, simple and efficient treatment concept
- Elegant, well thought-out instruments

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34

4628/1.000

Stiftset TMC Extrusion
Set "TMC Extrusion Pin"

97502L15.000.120	1	
97503.000.120	1	
97505L25.000.120	1	
97500.000.032	10	
97500.000.064	10	
97501.000.032	10	
97501.000.048	10	
97501.000.064	10	

Je 1 Stift/Linse/Stange, 5 x 10 Spanngummis
1 pin/lenticular element/bar each, 5 x 10 intraoral elastics

4629.000


Instrumentenset TMC Extrusion
Instrument set TMC Extrusion

S6881.314.012	1	
953M.314.014	1	
953AM.314.014	1	
8390.314.014	1	
H254E.314.012	1	
943CHZ.204.080	1	

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97505L25.000. ...

120

TMC Extrusionsstange/Steg aus glasfaserverstärktem Composite
TissueMaster Concept
TMC Extrusion bar made of glass fibre reinforced composite
TissueMaster Concept

97501

		10	10	10	10
Größe · Size		032	048	064	095

97501.000. ...

032

048

064

095

Intraorale Gummiringe, strong
Silikon, latexhaltig
Intraoral Elastics, strong pull
Silicone, containing latex

97500

		10	10	10
Größe · Size		032	064	095

97500.000. ...

032

064

095

Intraorale Gummiringe, medium
Silikon, latexhaltig
Intraoral Elastics, medium pull
Silicone, containing latex

Dr. Stefan Neumeyer

MaxilloPrep

Bone spreading and bone condensing

Developed in cooperation with Dr. Stefan Neumeyer, the **MaxilloPrep Spread and Condense screws** made of pure titanium guarantee effective, yet conservative bone spreading and condensing as part of treatments involving implants. These screws mobilize the bone during expansion and facilitate the creation of the bone bed by making use of the residual elasticity of the bone.

- The Spread and Condense screws are available in 5 sizes. They are screwed in and out in ascending order, until the required target size is reached. The sizes 4 and 5 also come with a triangular cross-section for flexible spreading of extremely hard bones.

- Produced according to up-to-date clinical know-how and the latest production techniques, these screws are particularly comfortable to use. Thanks to their double thread, these screws can be screwed in and out in no time at all. Asymmetrical thread flanks with a smaller angle in inserting direction facilitate the insertion of the screws. Penetration into the bone may also be assisted by percussion, i.e. by gently hammering the screw into the bone. To further facilitate the work of the operator, an optional mechanical wrench and a ratchet for use with the manual wrench are also available.

- The screws come without o-rings, which makes them particularly easy to clean. A clearly audible clicking sound confirms that the screw has securely locked into the hand wrench.

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

Bone Splitting, Bone Spreading with round and triangular screws

38

4596.000


Set MaxilloPrep Spread-Condense
Set MaxilloPrep Spread-Condense

	📄	
558L15.000.1	1	
558L15.000.2	1	
558L15.000.3	1	
558L15.000.4	1	
558L15.000.5	1	
559L15.000.4	1	
559L15.000.5	1	
186A.204.1	1	
210L16.205.013	1	
210L16.205.018	1	
538.000.	1	
538.204.	1	
534.000.	1	
537.000.	1	

Je 1 Spread-Condense-Schraube, Instrumentarium, Sterilcontainer, Inserttray
1 Spread-Condense screw, instrument set, sterilisation container and insert tray each

558 L 15

	📄	1	1	1	1	1
Größe · Size		1	2	3	4	5
L	mm	15	15	15	15	15

558L15.000. ...

	1	2	3	4	5
--	---	---	---	---	---

Gebrauchsmuster, Patente/ Utility model, patents
DE 10 2009 060 656*
* angemeldet/ * pending

MaxilloPrep Spread-Condense Schraube

Gr.1: 1,6-2,4 mm

Gr.2: 2,2-3,0 mm

Gr.3: 2,8-3,6 mm

Gr.4: 3,4-4,2 mm

Gr.5: 4,0-4,8 mm

Reintitan

MaxilloPrep Spread-Condense screw

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


Dr. Ernst Fuchs Schaller

Angle Modulation System

Horizontal widening of the alveolar ridge thanks to flexible mobilization from the inside

The Angle Modulation System developed in close cooperation with Dr. Ernst Fuchs Schaller allows a minimally invasive horizontal widening of the alveolar ridge even in those cases where there is not enough bone available to successfully insert an implant. The dental surgeon is frequently faced with an initial knife-ridge situation, especially in the lower jaw, in which case the widening of the alveolar ridge is an indispensable precondition.

The new surgical angle modulation method is based on the natural formation of new bone. The basic idea: The cortical lamella is flexibly mobilized from the inside and then successively expanded and erected towards the outside. An essential advantage of this non-traumatic method is that the jaw does not have to be folded up, i.e. the creation of a flap can be avoided. Consequently, the periosteum and the mucosa do not have to be detached and the "biologically active container" remains intact - to promote better healing without complications.

The essential advantage: Expansion of the alveolar ridge and insertion of the implant in just one session.

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
Diamantiert, konisch
Rostfreier Stahl
Diamond coated, tapered
Stainless steel

46 **4603.000**

Set Angle Modulation System A4 Schall
Set Angle Modulation System A4 with sonic tips

M110.000.016	1	
M110.000.028	1	
M110.000.034	1	
M110F.000.028	1	
M000FC.000.028	1	
574.000.	1	
H254E.314.012	1	
859.314.010	1	
859.204.010	1	
210L16.205.013	2	
SFS100.000.	1	
SFS101.000.	1	
SFS110.000.	1	

454.000.	1	
SF1979.000.	1	
566.000.	1	
186A.204.1	1	
210L16.205.018	1	
538.000.	1	
538.204.	1	
558L15.000.1	1	
558L15.000.2	1	
558L15.000.3	1	
558L15.000.4	1	
558L15.000.5	1	
559L15.000.4	1	
559L15.000.5	1	
TRC.205.055	1	
TRT.205.055	1	
9126.204.042	1	

Inkl. Inserttrays, Sterilcontainer A4
Including insert trays, sterilisation container A4

4604.000

Set Angle Modulation System A4
Set Angle Modulation System A4

M110.000.016	1	
M110.000.028	1	
M110.000.034	1	
M110F.000.028	1	
M000FC.000.028	1	
574.000.	1	
H254E.314.012	1	
859.314.010	1	
859.204.010	1	
210L16.205.013	2	

454.000.	1	
186A.204.1	1	
210L16.205.018	1	
538.000.	1	
538.204.	1	
558L15.000.1	1	
558L15.000.2	1	
558L15.000.3	1	
558L15.000.4	1	
558L15.000.5	1	
559L15.000.4	1	
559L15.000.5	1	
TRC.205.055	1	
TRT.205.055	1	
9126.204.042	1	

Inkl. Inserttrays, Sterilcontainer A4
Including insert trays, sterilisation container A4

48 **4605.000**

Set Angle Modulation System A6 Schall
Set Angle Modulation System A6 with sonic tips

M110.000.016	1	
M110.000.028	1	
M110.000.034	1	
M110F.000.028	1	
M000FC.000.028	1	
574.000.	1	

H254E.314.012	1	
859.314.010	1	
859.204.010	1	
210L16.205.013	1	
SFS100.000.	1	
SFS101.000.	1	
SFS110.000.	1	
SF1979.000.	1	
566.000.	1	

Inkl. Inserttrays, Sterilcontainer A6
Including insert trays, sterilisation container A6

4606.000

Set Angle Modulation System A6
Set Angle Modulation System A6

M110.000.016	1	
M110.000.028	1	
M110.000.034	1	
M110F.000.028	1	
M000FC.000.028	1	

Inkl. Inserttrays, Sterilcontainer A6
Including insert trays, sterilisation container A6

574.000.	1	
H254E.314.012	1	
859.314.010	1	
859.204.010	1	
210L16.205.013	1	

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Insert tray Angle Modulators, PP with silver ions, antibacterial (without instruments)

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Insert tray module frame, PP with silver ions, antibacterial (without instruments)

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Insert tray handles and manual wrenches, PP with silver ions, antibacterial (without instruments)

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Insert tray sonic tips, PP with silver ions, antibacterial (without instruments)

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Prof. Dr. Klaus-Ulrich Benner

Pipe Augmentation System

Augmentation of the alveolar ridge by depositing bone forming augmentation material underneath the periost

Developed by Prof. Dr. Klaus Ulrich Benner and Dr. Dr. Karl-Heinz Heuckmann, the pipe technique is a surgical method for increasing the height and the width of the alveolar ridge by depositing bone forming augmentation material underneath the periost. Each step of the process is carried out with the perfectly designed instruments of the Pipe Augmentation System, right down to the filling of the pipe with augmentation material. The material hardens upon contact with body fluid.

The pipe technique is particularly suitable for increasing the height and the width of atrophied alveolar ridges as a pre-prosthetic measure and prior to the insertion of an implant. As a result, sufficient bone is formed to provide an adequate bone structure.

Advantages:

- Minimally invasive augmentation method
- The bone forming material is moulded immediately after insertion into the pipe
- The method requires little time
- Merely a tiny insertion across the crest is required
- Minimal trauma during the intervention

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

4626.000

Set Pipe Augmentation System A6
Set Pipe Augmentation System A6

601.000.	1		
602.000.	1		
574.000.	1		
141P.106.027	1		
603.000.	1		
604.000.	1		
605.000.	1		
606.000.	1		
568.000.	1		

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MaxilloPrep Bone

Bone grafting and congruent augmentation

The MaxilloPrep Bone centring burs TRC and trepan burs TRT are indicated for the safe preparation of bone cylinders. The inventor, Dr. Stefan Neumeyer from Eschlkam, Germany, also uses a combination of these instruments to carry out a systematic and extensive bone augmentation, because the three diameters 045, 055 and 065 are ideally coordinated, allowing the fixation of bone cylinders of matching shapes. Dr. Neumeyer's basic idea is an integral evaluation of both the hard bone substance and soft tissue, in order to guarantee an optimum aesthetical clinical result. The line of MaxilloPrep instruments offers optimally adapted instruments and accessories to suit this comprehensive range of treatment.

- *The TRC centring bur creates a circular guide groove. The centring tip of the bur enables exact and safe positioning on the bone substance.*

- *The TRT trepan bur featuring a state-of-the-art blade configuration, arranged in pairs, is used for the preparation of the bone cylinder. The large perforations in the working part allow both unobstructed view of the operative site and easy removal of the bone cylinders.*

The innovative bone chip extractor 9126 as suggested by Dr. Stefan Neumeyer is used for precise, yet gentle removal of bone chips.

- *After exposing the bone substance by means of a very thin cut, the bone chip extractor is axially positioned on the bone. The self-centring tip penetrates the bone substance and the funnel-shaped bone chip extractor fills with moist bone chips.*

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*Preparation of bone cylinders,
creating a circular guide groove*

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*Creation of the bone bed using centring
burs with increasing or decreasing diameter*

**Minimal-invasive
vertikale Augmentation**

Minimally invasive vertical augmentation

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*Removal of bone chips from
spongy bone substance*

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




Dr. Martin Dürholt

MaxilloPrep Bonefix

Osteosynthesis screws for fixing bone grafts

Prior to carrying out a restoration with implants, the implantologist is frequently faced with a situation where there is not enough bone available. Developed in cooperation with Dr. Martin Dürholt, the MaxilloPrep Bonefix osteosynthesis screws allow exact placing and fixing of the bone graft at the recipient site.

The system comprises high-grade screws made of pure titanium in size 1.5 mm. They are available in lengths of 6, 9 and 12 mm. These osteosynthesis screws are provided with a continuous fast-cutting thread which ends directly at the flat screw head. Secure hold in the handle adapter is guaranteed by the hexagonal reception with a standard width across flats of 1.2 mm. The attachable handle is particularly light in weight and perfectly ergonomic.

Excellent vision during re-entry thanks to the blue anodised surface of the screw.

Advantages

- Small number of components
- "Minitray" suitable for sterilisation
- Protected removal of screws without contamination
- Optimum recess module and high-quality sterilisation container
- Can be used in combination with existing burs
- Compatible with conventional implant screw drivers
- Suitable for all established bone augmentation techniques

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
577

MaxilloPrep Bonefix Handgriffadapter
Rostfreier Stahl
Maxillo Prep Bonefix handle adapter
Stainless steel

582


Abmessungen · Dimensions mm 46 x 14 x 17

Inserttray (Minitray) MaxilloPrep Bonefix
PP mit Silberionen, antibakteriell (ohne Instrumentarium)
Insert tray (mini tray) MaxilloPrep Bonefix
PP with silver ions, antibacterial (without instruments)

Empfohlenes Instrumentarium:
Pilotbohrer 210L16.205.013, 018,
Knochenfräser rund
H141A.104/205.027.

Recommended instruments:
pilot drill 210L16.205.013, 018,
round bone cutter
H141A.104/205.027

210L16.205.013

210L16.205.018

⊖_{opt.} 800 - 1 000 min⁻¹ / rpm

⊖_{max.} 6 000 min⁻¹ / rpm

H141A.205.027

H141A.104.027

⊖_{opt.} 6 000 min⁻¹ / rpm

⊖_{max.} 100 000 min⁻¹ / rpm

4609.000

Set MaxilloPrep Bonefix Osteosyntheseschrauben
Set MaxilloPrep Bonefix osteosynthesis screws

583L6.000.015	5		
583L9.000.015	5		
583L12.000.015	5		
577.000.	1		
574.000.	1		
582.000.	1		

Inkl. Inserttray (Minitray)
Including insert tray (mini tray)

4610.000

Set Sterilcontainer A8 (90 x 90 x 55 mm) und Inserttray mit Mulden (ohne Instrumentarium)
Set sterilisation container A8 (90 x 90 x 55 mm) and insert tray with recesses (without instruments)

571.000.	1		Inserttray mit Mulden Insert tray with recesses
556.000.	1		Sterilcontainer A8 Sterilisation container A8

Universelle Pilotbohrer Implantologie

Universal pilot burs for use in Implantology

The universal pilot burs are ideally suited for the initial preparation of the axis and depth of an implant site. If an implant with large diameter is to be placed, it is recommended to enlarge the perforation gradually.

The pilot burs with particularly small diameter can be used for palpating the implant position, perforation of bone cylinders as well as for creating a continuous perforation line during bone spreading.

Advantages:


- Large chip spaces for good chip removal
- Pyramid-shaped, special instrument tip for easy penetration
- Effective cutting
- Lasered depth markings at intervals of 2 mm, starting at 8 mm from the instrument tip
- Size and instrument length lasered on to the shank for easy identification

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





8, 10, 12, 14 (16, 18) mm bei Gr. 020-028


Pilot bur for implantology, stainless steel

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Tiefenmarkierung = 8, 10, 12, 14 (16, 18) mm

Pilot bur for implantology, made of ceramics

Depth marking = 8, 10, 12, 14 (16, 18) mm


Stainless steel bur blocks

A range of functional bur blocks is available for cleaning, disinfection, storage and sterilization of dental instruments.

The blocks are provided with silicone plugs to guarantee secure hold of the instruments and to prevent them from falling out.

All stainless steel bur blocks are suitable for cleaning and disinfecting in the instrument or ultrasonic bath or in the thermo disinfecter. Komet had the reprocessing of our instruments validated by an external institute. With Komet bur blocks and reprocessing instructions, our customers are always on the safe side.

- *Different types of shank can be combined. The bur block fits turbine, hand-piece and contra-angle instruments.*
- *Equally suitable for organizing preparation sequences, e.g. rotary instruments to perform a sinus floor elevation.*
- *9890L4 / Height 40 mm, suitable for example for bone cutters H254E/LE.314, Trepan burs 227A/B/C.204*
- *9890L5 / Height 50 mm, suitable for example for bone cutters H141A, H162A, H162ST, H166A.104, Pilot drills for implantology 210L16/L19/L20.204 and .205*
- *9890L7 / Height 70 mm, suitable for instruments with shank 105, long hand-piece*

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Bur block made of stainless steel with 16 blue silicone plugs as universal instrument holders, for FG and RA instruments with a maximum length of 33 mm

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Universal bur block made of stainless steel with 6 blue silicone plugs as universal instrument holders, different types of shanks can be combined, for a maximal instrument length of 37 mm

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Universal bur block made of stainless steel with 6 blue silicone plugs as universal instrument holders, different types of shanks can be combined, for a maximal instrument length of 47 mm

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Universal bur block made of stainless steel with 6 blue silicone plugs as universal instrument holders, different types of shanks can be combined, for a maximal instrument length of 67 mm

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Bur block made of stainless steel with 12 blue silicone plugs for FG and RA instruments with a maximum length of 58 mm

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Bur block made of stainless steel with 24 blue silicone plugs for FG and RA instruments with a maximum length of 58 mm

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Bur block made of stainless steel with 40 blue silicone plugs as universal instrument holders, for FG and RA instruments with a maximum length of 58 mm

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Bur block made of stainless steel with 6 blue silicone plugs for FG and RA instruments and 3 plugs for sonic tips, for a maximum length of 58 mm

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Insert trays and Sterilization container

Systematic insert trays and safe sterilization containers

Many systems of the Komet® dental range comprise handy insert trays for the reception of instruments and system components. All components are clearly and logically arranged in the correct sequence ready for use during treatment sessions. Both insert tray and instruments are completely integrated in a high-quality sterilization container.

Insert tray

- The insert tray is made of a special high-performance plastic.
- Antimicrobial effect thanks to the continuous release of silver ions
- High purity grade during idle times
- All system components are clearly arranged.

Sterilization container

- Safe sterilization. Does not require sealing in foil
- No seals, no maintenance required. Lasts through more than 2,000 sterilization cycles
- PTFE filter membrane with a particularly small nominal pore size of 0.2 µm – retains even airborne germs and aerosols
- Easy filter change after 150 cycles

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Osteosynthesis screws made of pure titanium

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The bone graft, perforated with the pilot drill 210L16.

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Countersinking of the screw heads - preparation of the perforations with the round bur H141A.104.027.

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The bone graft is screwed into the cortical bone of the recipient site with one of the screws of the MaxilloPrep Bonefix system.

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Minimally invasive widening of the alveolar ridge

Typische Knife-ridge Situation

Typical initial situation: Knife-ridge

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Incision and exposure, instrument is slightly tilted in lingual direction, parallel to the defect

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Mobilization of the bone: Deepening the inner side of the buccal lamella and creating vertical incisions at the mesial and distal end of the cortical bone.

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Slow, successive spreading whilst continually checking the flexibility with the flat modulator.

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Erection of the buccal cortical lamella, insertion of the implant during the same treatment session.

Komet Dental

Gebr. Brasseler GmbH & Co. KG

Trophagener Weg 25 · 32657 Lemgo

Postfach 160 · 32631 Lemgo · Germany

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Telefon +49 (0) 5261 701-700

Telefax +49 (0) 5261 701-289

info@kometdental.de

www.kometdental.de

Export:

Telefon +49 (0) 5261 701-0

Telefax +49 (0) 5261 701-329

export@kometdental.de

www.kometdental.de

Komet Austria Handelsagentur GmbH

Innsbrucker Bundesstraße 75

5020 Salzburg · Austria

Telefon +43 (0) 662 829-434

Telefax +43 (0) 662 829-435

info@kometdental.at

www.kometdental.at