MAJOR PLUS / COMP

**WFA**

WIDE FUNCTIONAL AREA

MOULD CHART

MAJOR

# MAJOR PLUS \* MAJOR PLUS COMP \*\*

# WFA WIDE FUNCTIONAL AREA

The new MAJOR PLUS WFA tooth follows the mounting philosophy according to Prof. Albert Gerber, the palatal cusp of the upper 5th, 6th and 7th tooth occludes in the respective fossa of the lower 5th, 6th and 7th tooth; in the 4th tooth this occlusion is reversed, the vestibular cusp of the lower tooth occludes in the fossa of the upper 4th tooth. The term MAJOR PLUS WFA (Wide Functional Area) indicates the type of occlusion done during the mastication, when not only one point of contact is affected but a large

functional area of around 2 mm.

## OCCUSAL PROFILES

The occlusal tables respect the concept mortar and pestle. The wide fossa allows a large spectrum of occlusal contacts without interferences. MAJOR PLUS WFA teeth have such an high precision to reduce to the minimum any possible selective grinding made by the operator.

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
## \* MAJOR PLUS WFA

Major Plus is a multi-layer acrylic tooth designed to comply with the higher functional and aesthetic standards of modern prosthetic technique.

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## MOULDS AND SHADES

The MAJOR PLUS WFA tooth is available in the universal A1-D4 shade guide. The need to simplify and speed up the clinical and technical procedures, maintaining high quality standard, suggested the selection of three posterior moulds only: SMALL, MEDIUM, LARGE, adaptable to any of the anterior MAJOR PLUS moulds.

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EN

SYNTHETIC POLYMER TEETH FOR DENTAL PROTHESIS

- Store in a clean, dry place away from light.
- To ensure optimum bonding to the denture base resin, grind teeth on the cervix prior to setting up, remove wax and degrease with a methyl methacrylate monomer.
- Direct the patient to adopt daily hygiene protocols, suggest him to avoid excessive brushing which could cause abrasive wear, and not to immerse in alkalis, acids or other substances that can damage acrylic materials.
- This product does not contain substances classified as hazardous conforming to Directive 67/548/EEC and subsequent amendments. Consequently, no Material Safety Data Sheet is available.
- Grinding of the teeth generates dust which is classified as "not otherwise classifiable". Do not inhale. Use appropriate protection. TLV/TWA applies.
- Dispose conforming to local law. Product and packing materials are not listed in the Hazardous Waste Inventory as per C.D. 91/689/EEC 1(4). Patient-contacted materials are at infection risk even for disposal purposes (EU Hazardous Waste Inventory No. 180103).
- Available in the following shade guides: A1-D4.

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MAJOR Prodotti Dentari S.p.A.